



Product Characteristics

The most speedy and wash disinfection
Touch-screen so that the program parameters can be displayed and controlled
Auto door with foot switch
Automatic visual glass door adopted for observing the operation inside the tank
Totally enclosed to avoid sanitizer odor leakage
Automatic adjuvant injection: only need to link the adjuvant container with the tank through a connector

Work with various kinds of sanitizer: glutaraldehyde, peracetic acid, phthalic dicarboxaldehyde, chlorine dioxide and acidic electrolyzed oxidizing water.
Includes leaking test, blocking test, wash, enzyme wash, rinse, disinfection, final rinse, alcohol drying and air drying
Data management and printing for disinfection information
Overall dimensions: 74x65x94
Weight 100Kg
Water Consumption/cycle: 25L
Dynamic drainage

AFG

Avenida do Infante, n 8
Edifício Executivo, 1F
9000 - 015 Funchal, Madeira

291 763 044
geral@afg.pt
www.afg.pt