



### Product Characteristics

Incubator for many applications, such as, plant propagation, humidity tests, simulation of tropical conditions, incubation of complete instruments, corrosion testing;  
Full visibility of samples thanks to the acrylic structure;  
Easy access to working chamber;  
Forced air circulation and electronic temperature control;  
Temperature range: from +5°C to 60°C;  
Capacity: 60 litres;  
Internal dimensions: 450 x 380 x 380mm;  
External dimensions: 600 x 390 x 390mm;  
Net Weight: 11,2Kg;  
Power: 230V, 50Hz, 350W;  
Optional: digital model (SI60D) and analogue model (SI60).

#### AFG

Avenida do Infante, n 8  
Edifício Executivo, 1F  
9000 - 015 Funchal, Madeira

291 763 044  
geral@afg.pt  
www.afg.pt