



Product Characteristics

Incubator for many applications, such as, plant propagation, humidity tests, simulation of tropical conditions, incubation of complete instruments, corrosion testing;
Full visibility of samples thanks to the acrylic structure;
Easy access to working chamber;
Forced air circulation and electronic temperature control;
Temperature range: from +5°C to 60°C;
Capacity: 60 litres;
Internal dimensions: 450 x 380 x 380mm;
External dimensions: 600 x 390 x 390mm;
Net Weight: 11,2Kg;
Power: 230V, 50Hz, 350W;
Optional: digital model (SI60D) and analogue model (SI60).

AFG

Avenida do Infante, n 8
Edifício Executivo, 1F
9000 - 015 Funchal, Madeira

291 763 044
geral@afg.pt
www.afg.pt