



Product Characteristics

- High flow rate accuracy is achieved: within $\pm 2\%$;
- Dedicated Mindray infusion sets be used delivering high accuracy over long periods;
- More than 9 hours of battery operating time with only 6 hours charging time;
- Various programming modes provide flexible options for clinical applications in different departments;
- TIVA mode, hands-on and hands-free Bolus setting makes an anesthesia infusion simple with precise delivery;
- All important parameters can be viewed from over 5 meters away, making hospital monitoring efficient and convenient;
- Recognize various syringe brands and sizes;
- Sizes: 5/10/20/30/50/60ml;
- Horizontal design allows the pump to be easily stacked and inserted into the Dock;
- Three levels of alarm for all security alarms (H, M and L) ,8 levels of adjustable alarm volume;
- Alarm light is large and intuitive,clearly visible over long distance to maxmize patient safety;
- Dimensions: 295mm (w) x 87mm (h) x 174 (d).

AFG

Avenida do Infante, n 8
Edifício Executivo, 1F
9000 - 015 Funchal, Madeira

291 763 044
geral@afg.pt
www.afg.pt